

Annexure-A

1	<p>Neonatal Nursing & CPR Manikin</p> <p>Infant anatomic characteristics, soft skin made of imported material, realistic touch feeling, flexible joints and capable of providing full range of nursing exercise.</p> <p>Features:</p> <ul style="list-style-type: none">• Pupil observation: Pupil states: normal, mydriasis and myosis etc.• Hair combing and face washing• Otic drops and irrigation• Oral cavity care• Airway management: realistic mouth, nose, tongue, gum, pharynx, larynx, esophagus, epiglottis, trachea and tracheal ring, can practice tracheal intubation sputum suction, oxygen inhaling• Oral and nasal tracheal intubation• Incision of tracheal nursing• Venipuncture, injection and transfusion in arm• Intramuscular injection in bilateral deltoid, bilateral vastus lateralis and bilateral buttock• Enema• Urethral catheterization• Holistic care: sponge washing, replacing clothes, cold and hot therapy• Limbs joints: bend, rotation, upper and lower movement• CPR operation training: should support many inhaling ways, e.g. mouth to mouth, mouth to nose, simple respirator to mouth; electronic monitoring of airway open, inhaling times, frequency and volume and compression times, frequency, site and depth; automatically judge the ratio of artificial respiration and external chest compression; show data in real time with English voice prompts during the entire process.
---	--



2	<p align="center">Child airway management Trainer with Stand</p> <p>Should include head mounted on sturdy stand, lungs, stomach, pump spray lubricant, instruction manual, and hard carry case.</p>
3	<p align="center">Complete Child Airway Management Kit</p> <p>The Complete Child Airway Management Kit including Smith's selection of manual resuscitator with mask, three 4.0 mm uncuffed ET tubes, size 2.5 larynx mask airway, and one each of Guedel, Berman and nasal airways; laryngoscope set that includes the handle and sizes 2, 3 and 4 Macintosh blades and sizes 0, 1, 2 and 3 Miller blades; stylette; and lube (6 oz). Manikin should include head mounted on sturdy stand, lungs, stomach, and pump spray lubricant.</p>
4	<p align="center">Total Adult Airway Management Kit</p> <p>Total Adult Airway Management Kit containing Smith's selection of one Nu-trake kit, manual resuscitator with mask, three 7.5 mm uncuffed ET tubes size, 5 larynx mask airway, one each of Guedel, Berman and nasal airways; laryngoscope set that includes the handle and sizes 2, 3, and 4 Macintosh blades and sizes 0, 1, 2 and 3 Miller blades; stylette; and lube (8 oz.) size: 12"x 10"x10". Sh. Wt. 5 lbs.</p>
5	<p align="center">Nursing Manikin Female</p> <p>Full-functional Elderly Nursing Manikin (Female):-</p> <p>Basic nursing functions:</p> <ul style="list-style-type: none"> • Hair and face washing • Pupil observation, Bilateral pupil states; normal, mydriasis • Eye and ear washing and administering • Removal and insertion of hearing aid • Oral cavity and artificial teeth care • Endotracheal intubation, automatically monitor whether the tracheal intubation is correct or not through sensor, support auscultation to detect intubation site. • Incision of trachea nursing • Sputum suction • Oxygen inhaling • Oral and nasal feeding • Gastrolavage • Venipuncture, injection, transfusion (arm) • Deltoid subcutaneous injection • Vastus lateralis injection • Thoracic cavity, abdominal cavity, liver, bone marrow and lumbar puncture • Breast nursing, mammary gland examination • Enema • Male/Female urthelial catheterization • Male/Female bladder irrigation



7

- Standard American Heart Association 2010 guideline for CPR and ECC, to comply with

CPR manikin features:

- Allow to perform airway open, artificial breathing and external cardiac compression.
- While performing artificial aspiration, judge inflation volume by chest rise
- External chest compression, Correct compression depth: correct compression site, prompt correct buzzing sound Excessive compression site, prompt alarming buzzing sound
- Operation frequency: >100 / min, with audio
- Battery powered.

Schlack